

t:slim X2™ Insulin Pump

t:lock™ Infusion Sets & Cartridge



AutoSoft™ 90




- 90° insertion angle
- Integrated inserter
- Teflon cannula

AutoSoft™ 30



- 30° insertion angle
- Integrated inserter
- Teflon cannula

TruSteel™ Set



- Manual straight in insertion
- Stainless steel/metal cannula

VariSoft™ Set



- Manual angled set
- Teflon cannula

t:lock Infusion Set NDSS Order Codes (10 per pack)

Tandem AutoSoft 90 Infusion Set t:lock 6mm/60cm	826
Tandem AutoSoft 90 Infusion Set t:lock 6mm/110cm	827
Tandem AutoSoft 90 Infusion Set t:lock 9mm/60cm	828
Tandem AutoSoft 90 Infusion Set t:lock 9mm/110cm	829
Tandem AutoSoft 30 Infusion Set t:lock 13mm/60cm	830
Tandem AutoSoft 30 Infusion Set t:lock 13mm/110cm	831
Tandem VariSoft Infusion Set t:lock 13mm/60cm	832
Tandem VariSoft Infusion Set t:lock 17mm/60cm	833

Tandem TruSteel Infusion Set t:lock 6mm/60cm	834
Tandem TruSteel Infusion Set t:lock 6mm/80cm (not currently available)	835
Tandem TruSteel Infusion Set t:lock 8mm/60cm	836
Tandem TruSteel Infusion Set t:lock 8mm/80cm (not currently available)	837

t:lock Cartridge

Tandem t:lock Cartridge	722
-------------------------	-----

Cannula Fill and Site Change Guidelines for Infusion Sets

	Fill	Site Change
AutoSoft 90	6mm: 0.3 units, 9mm: 0.5 units	2-3 days
AutoSoft 30	0.7 units	2-3 days
TruSteel	0 units/no fill	1-2 days
VariSoft	0.7 units	2-3 days

For more information on the t:slim X2 Insulin Pump, please contact the AMSL Diabetes Customer Care Team on **1300 851 056** or at diabetes@amsl.com.au

amsl diabetes.com.au    

Always read the label and use only as directed. Read the warnings available on www.amsl diabetes.com.au/resources before purchasing. Consult your healthcare professional to see which products are right for you. ARTG 116234, ARTG 304681. PR-100-50 January 2019

